

Public Health England Microbiology Services Porton Down Salisbury Wiltshire SP4 OJG

## **Certificate of Containment Testing**

## Containment Testing of Caps for Rotor S-4-104 with DWP-Buckets in the Eppendorf 5810/R Bench Top Centrifuge

Report No. 111/13 A

Report Prepared For: Eppendorf AG, Hamburg, Germany

**Issue Date:** 

10<sup>th</sup> April 2014

## **Test Summary**

Caps for rotor S-4-104 with DWP-Buckets were containment tested in the Eppendorf 5810/R bench top centrifuge, using Annex AA of IEC 61010-2-020:2006 (2<sup>nd</sup> Ed.). The sealed rotor was shown to contain a spill within the centrifuge.

**Report Written By** 

Name: Miss Anna Moy

Title: Biosafety Scientist

**Report Authorised By** 

Name: Mrs Sara Speight

Title: Senior Biosafety Scientist

ps.M.

Please be aware that the use of the Royal Coat of Arms is highly restricted and cannot be copied. Please do not put the PHE logo on your website or use our name to endorse your products. Any reference to PHE needs to be approved by us before it can be used.