


Certificate of Purity - BIOPUR®

This package contains a high-quality consumable manufactured under the Biopur Eppendorf Purity Standard. The Eppendorf Biopur consumables are produced in a class 6 (according to VDI 2083) and a class 100,000 (according to U.S. Fed. Stand. 209 D) clean room environment.

For this product Eppendorf certifies the following:

- > Sterile
- > Pyrogen-free
- > Human DNA-free
- > Bacterial DNA-free
- > DNase-free
- > RNase-free
- > PCR inhibitor-free
- > ATP-free



These parameters are continuously monitored by an independent certified laboratory. Eppendorf guarantees the conformity within the following limits:

Sterility	in accordance with USP, Ph. Eur. 2.6.12
Pyrogens (Endotoxins)	< 0.001 EU/ml (tested according to Ph. Eur. 2.6.14 (LAL test))
Human DNA	< 2 pg
Bacterial DNA	< 50 fg
DNase	not detectable (LOD: 1.0×10^{-6} Kunitz units)
RNase	not detectable (LOD: 1.0×10^{-9} Kunitz units)
PCR inhibition	less than 10 targets amplifiable
ATP	< 5.5×10^{-12} mg

Lot-specific certificates can be downloaded from the Internet at www.eppendorf.com/certificates.

March, 2016

Page 1

Daniel Voss
Product Line Manager
Consumables

Dr. B. Schreiber
Vice President Quality Management &
Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf®, the Eppendorf logo and Biopur® are registered trademarks of Eppendorf AG, Germany. U.S. Design Patents are listed on www.eppendorf.com/ip. All rights reserved, incl. graphics and pictures. Copyright 2016 © by Eppendorf AG.