Health Protection Agency Microbiology Services Porton Down Salisbury Wiltshire SP4 0JG



## Certificate of Containment Testing

Containment Testing of Rotor S-4-104 with Round Buckets (5820 741.007-00) in the Eppendorf 5810/R Bench Top Centrifuge

**Report No. 196-12 A** 

Report Prepared For: Eppendorf AG, Hamburg, Germany

**Issue Date:** 12<sup>th</sup> September 2012

## **Test Summary**

Rotor S-4-104 with Round Buckets (5820 741.007-00) was containment tested in the Eppendorf 5810/R bench top centrifuge, using Annex AA of IEC 1010-2-20. The sealed rotor was shown to contain a spill within the centrifuge

**Report Written By** 

Anna

/

Name: Miss Anna Moy

Title: Biosafety Scientist

Name: Mrs Sara Speight

Report Authorised By

**Title: Senior Biosafety Scientist**